

gibco



Gibco cell culture reagents

Designed to deliver reproducibility and performance
for results you can count on every day

ThermoFisher
SCIENTIFIC



Gibco media—over 50 years of cell culture expertise

Gibco™ media, sera, supplements, cells, and cell culture reagents are designed to deliver reproducibility and performance for results you count on.

Regardless of whether you are performing cell culture in a research lab or production facility, need a special or defined formulation, or are growing cell lines, primary cells, or stem cells, Gibco products offer a reliable solution.

Contents



Cell culture reagents	
Media	4
Media: stem cells	6
Supplements	7
Sera	8
Dissociation reagents	10
Balanced salts	11
Growth factors	12
Antibiotics	13
Cryopreservation	14
Matrices	16
<hr/>	
Cell culture	
Custom media	17
Equipment and instrumentation	18
<hr/>	
Gibco Education	19
<hr/>	
Culture then transfect	20
<hr/>	
Culture then express	21
<hr/>	

Cell culture reagents—media

Time-tested and trusted, Gibco media include products designed to support the growth and maintenance of a variety of mammalian cells and cell lines

We've developed ready-to-use media products as well as powdered and concentrated liquid formulations to fit your experimental setup and budget. Classical media products include DMEM, RPMI 1640, MEM, IMDM, and Opti-MEM™ I media.

Classical media product	Summary	Web link
DMEM	Growth of a spectrum of mammalian cell lines	thermofisher.com/dmem
RPMI 1640	Enriched media with extensive applications for mammalian cells	thermofisher.com/rpmi
MEM	Patterned after Eagle's Media, well suited for mammalian cells	thermofisher.com/mem
IMDM	Highly enriched media well suited for rapidly proliferating, high-density cell cultures	thermofisher.com/imdm
Opti-MEM I medium	Ideal for use during cationic lipid transfections	thermofisher.com/optimem



In addition to these classical media, we have developed products to upgrade your cell culture:

Upgrade basal media to Advanced media, upgrade L-glutamine to GlutaMAX media, and upgrade PBS/DMEM to FluoroBrite DMEM.

Culture cells with Advanced media

Consistent cell growth, less serum, lower cost

Enriched with normal-serum constituents, Gibco™ Advanced™ media can be used in combination with your fetal bovine serum (FBS), and feature:

- Formulations that support common cell lines (e.g., Advanced DMEM, DMEM/F-12, MEM, RPMI 1640)
- Fewer lot-to-lot changes of serum, which means less variability
- Reduced cost due to lower serum usage and fewer lots tested
- Equivalent cell growth/no change in morphology of common cell lines (Figure 1)




Figure 1. Growth of cell lines in Advanced media. Growth comparison of multiple cell lines in standard RPMI 1640 supplemented with 5% FBS to Advanced RPMI 1640 supplemented with 0.5%, 2%, and 5% FBS.

Select your Advanced media at thermofisher.com/advanced

Culture cells with GlutaMAX media

More stable than L-glutamine—keeping your cells healthier for longer

L-glutamine can degrade in culture media to form toxic ammonia, which decreases cell viability. Gibco™ GlutaMAX™ Supplement is a stabilized form of L-glutamine, and features:

- Superior stability—helps prevent degradation during storage or in the presence of cells, which happens with L-glutamine (Figure 2)
- Flexibility—available as a stand-alone supplement or as a component in Gibco™ basal media products




Figure 2. Stability of Gibco™ CTS™ GlutaMAX™-I Supplement vs. L-glutamine in DMEM. DMEM was supplemented with GlutaMAX-I Supplement or L-glutamine, aliquotted into vials, and stored at 37°C. Samples were taken daily and frozen at -20°C. Levels of GlutaMAX-I Supplement and L-glutamine were determined by HPLC.

Select your GlutaMAX media or supplement at thermofisher.com/glutamax



Image cells with FluoroBrite DMEM

Reduces background fluorescence in cell imaging while preserving cell viability

Gibco™ FluoroBrite™ DMEM is a culture medium designed for imaging cells, and features:

- Background autofluorescence 90% lower than the phenol red-free DMEM equivalent (Figure 3)
- Improved signal-to-noise ratio compared to the phenol red-free DMEM equivalent




Figure 3. Fluorescence of PBS, FluoroBrite DMEM, and phenol red-free DMEM at 509 nm (excitation at 480 nm).

Select FluoroBrite DMEM at thermofisher.com/fluorobrite

Cell culture reagents—media: stem cells

Manage your pluripotent stem cell (PSC) culture schedule and eliminate daily feeding

Gibco™ Essential 8™ Flex Medium is formulated to extend the activity of key heat-sensitive components found in PSC medium, including FGF2, to enable a truly weekend-free culture feeding schedule.

- **Flexible feeding schedule**—maintain pluripotency over a full 2-day period without feeding cells

- **Easy to transition**—same setup, applications, and protocol as the original Essential 8 Medium formulation, just without the daily feeding
- **Proven**—based on the original Essential 8 Medium formulation, maintains pluripotency and normal karyotypes in long-term culture

The old way: feed your cells every day, 7 days a week



Weekend-free feeding is possible—compare your PSC feeding schedule to Essential 8 Flex Medium.

The Flex way: eliminate daily feeding



Learn more at thermofisher.com/essential8flex

Cell culture reagents—supplements

Gibco media supplements help ensure reliability and consistency of your cell culture research

Advantages of using media supplements include customizing the growth conditions of your cells, improving cell viability and growth, and keeping cells healthier longer.

Gibco media supplements include:

2-Mercaptoethanol	Insulin	Amino acids
N-2 Supplement	B-27 Supplement	Primary cell supplements
BSA	Sodium pyruvate	Glutamine or GlutaMAX media
Transferrin	G-5 Supplement	

Select your media supplement at [thermofisher.com/mediasupplements](https://www.thermofisher.com/mediasupplements)

Culture neural cells with Gibco B-27 Supplement

The most cited neuronal cell culture supplement

- Gibco™ B-27™ Supplements support the culture of primary neurons in low or high densities in both short- and long-term cultures
- Several additional formulations have been developed for specific customer applications, including stem cell differentiation and enhanced electrophysiology experiments
- For over 20 years, B-27 products have been manufactured under cGMP guidelines to help ensure reliability, traceability, and lot-to-lot consistency to minimize variability

Choose your B-27 Supplement at [thermofisher.com/b27](https://www.thermofisher.com/b27)



Cell culture reagents—sera

Gibco sera have earned the trust of researchers around the world due to consistent quality, superior confidence, and award-winning support

Gibco sera help meet your research needs and budget requirements, offering the best value for basic cell culture, specialty research, and specific assays. We are dedicated to providing you with the right fetal bovine serum (FBS) for your specific cell culture need and your lab budget.

Select your serum at thermofisher.com/fbs

Origin and traceability reassurance like no other

We continue to take the necessary steps to be the unrivaled quality leader in FBS

Over the past few years, we've provided researchers more transparency into the collection and manufacturing process of Gibco FBS. We consider it imperative that our customers know exactly what they're getting and have confidence that our products will enable the results they demand.

International Serum Industry Association (ISIA)
traceability certification

Effective February 12, 2014, ISIA traceability certification was awarded to Thermo Fisher Scientific

This certification provides researchers with peace of mind and the confidence that Gibco sera are manufactured under the highest traceability standards, offering quality and performance for your research.

Gibco FBS fingerprinting, powered by Oritain

We have developed a technology in partnership with Oritain that helps ensure the quality of sera obtained from low-viral risk regions

Gibco FBS fingerprinting powered by Oritain measures the naturally occurring elements in the environment that are absorbed by plants, animals, and soils, and then compiles those measurements into a "fingerprint." Using this fingerprint, we can link serum to its specific geographical location—providing a robust means to verify if a product is genuine.

Gibco sera are the first to be offered with this technology, giving scientists reassurance regarding the origin of each serum.

Find out more about ISIA certification and FBS fingerprinting at thermofisher.com/trustfbs



The right sera for all your cell culture needs

Gibco serum category*	Standard	Performance
	Sera for cell culture with robust cell lines—excellent value for basic research	Low-endotoxin sera for general cell culture with common cell lines
Recommended products	FBS Qualified, USDA-approved origins, South America, Canada	FBS Qualified, US
Endotoxin specification/standard	Typically ≤50 EU/mL	≤10 EU/mL
Quality and performance testing (including standard tests: growth, cloning, plating)	Standard testing	Standard testing, plus exclusive BVDV screening of raw material (FBS)
Popular catalog numbers/standard	10437028 FBS Qualified, USDA-approved origins 10270106 FBS Qualified, South America	26140079 FBS Qualified, US 16140071 FBS Qualified, heat inactivated, US

Gibco serum category*	Performance Plus	Secure**	Specialty
	Lowest-endotoxin and most highly characterized sera; good for broad range of cell types, especially sensitive cell lines	Sera sourced from BSE-negligible regions for preclinical, industrial, and academic research applications requiring low risk	Sera qualified for specialty research and specific assays, including stem cell research, immunoassays, antibodies, and others
Recommended products	FBS Certified, US	FBS Qualified, Australia FBS Qualified, New Zealand DBS, New Zealand	FBS and other sera for specialty research and assays
Endotoxin specification/standard	≤5 EU/mL	≤10 EU/mL for FBS, donor bovine sera (DBS)	Per Certificate of Analysis
Quality and performance testing (including standard tests: growth, cloning, plating)	Standard testing, plus analytical tests for hormone and biochemical profiles and exclusive BVDV screening of raw material (FBS)	Standard testing, plus exclusive BVDV screening of raw material (FBS) or donor animals (DBS)	Standard testing, plus exclusive BVDV screening of raw materials (FBS) or donor animals (DBS) prior to final manufacturing
Popular catalog numbers/standard	16000044 FBS Certified, US 10082147 FBS Certified, heat inactivated, US	10099141 FBS Qualified, Australia 10100147 FBS Qualified, heat inactivated, Australia 10091148 FBS Qualified, New Zealand	16141079 FBS ES Cell Qualified, US 12676029 FBS, Charcoal Stripped, USDA-approved origins

* All products may not be available in all regions due to importation regulations. Contact your local sales representative regarding product availability in your country.

** All Secure sera can be gamma irradiated upon request to comply with EU and US regulations and guidelines.

Cell culture reagents—dissociation reagents

Gibco cell dissociation products are ideal for use with tissues and cell monolayers. Trypsin and TrypLE reagents come in a wide variety of formats to meet the diverse needs of researchers performing adherent cell culture:

	Fast, general purpose	Gentle and convenient for research use	Gentle and convenient for bioproduction/industrial applications
	Trypsin	TrypLE™ Express (1X)	TrypLE™ Select (1X)
Source origin	Animal origin (porcine)	Animal origin-free (AOF)	Animal origin-free (dedicated AOF machinery)
Storage temperature	Frozen (-5 to -20° C)	Ready to use (room temp. stable)	Ready to use (room temp. stable)
Higher cell viability		Yes	Yes
Inactivation method	Requires trypsin inhibitors or neutralizers	Inhibition by dilution (no inhibitors needed)	Inhibition by dilution (no inhibitors needed)

Detach cells with TrypLE Express reagent

Cell detachment that is gentle on cell surface proteins

Gibco™ TrypLE™ reagents are highly purified, recombinant cell-dissociation enzymes that replace porcine trypsin. These reagents are ideal for dissociating attachment-dependent cell lines in both serum-containing and serum-free conditions. They can directly substitute for trypsin without protocol changes, and are:

- **Stable at room temperature**—no need to thaw
- **Gentle on cells**—help protect your cells' surface proteins (Figure 4)
- **Animal origin-free**—important if you need a product without animal-derived components

Select your TrypLE reagents at thermofisher.com/tryple




Figure 4. Comparison of cell dissociation reagents. Jurkat cells were used to demonstrate the effects of cell dissociation reagents on detection of the cell surface epitope CD2. **(A)** Time course of treatment with 1X Gibco™ TrypLE™ Express reagent. There is no reduction in the CD2 fluorescence signal. **(B)** Time course of treatment with 0.05% trypsin-EDTA. There is a time-dependent reduction in CD2 detection.

Discover the full range of dissociation reagents at thermofisher.com/celldissociation

Cell culture reagents—balanced salts

All Gibco balanced salt solutions are manufactured in state-of-the-art cGMP, ISO-certified facilities, helping ensure the highest quality and consistency for reproducible results

Balanced salt solutions can provide an environment that maintains the structural and physiological integrity of cells *in vitro*.

To ensure consistency, Gibco balanced salt solutions are tested for osmolality, pH, stability, the absence of bacterial and fungal contamination, and endotoxin. All powdered balanced salt solutions are produced without sodium bicarbonate to increase stability.



Our most popular balanced salts are below:

DPBS

DPBS, Dulbecco's Phosphate Buffered Salines

HBSS

HBSS, Hanks' Balanced Salt Solutions

PBS

PBS, Phosphate-Buffered Salines

EBSS

EBSS, Earle's Balanced Salt Solutions


Select your balanced salts at [thermofisher.com/balancedsalts](https://www.thermofisher.com/balancedsalts)

Cell culture reagents—growth factors

Select pure, high-quality growth factors to help you achieve consistent cell culture

Growth factors are the unique cell signaling molecules that play important roles in cell proliferation and development. Gibco growth factors are designed to give you:

- **High biological activity**—more results with less protein
- **High purity**—minimized interference from proteins or other contaminants
- **Proven compatibility**—Gibco proteins bioassayed with Gibco media



Some popular growth factors are below:

BDNF

Brain-Derived Neurotrophic Factor (BDNF) is used in the maintenance of motor neurons and the differentiation of both cholinergic and dopaminergic neurons.

bFGF

bFGF (also known as fibroblast growth factor-basic, FGF-basic, FGF-β, FGF2, or heparin binding growth factor) is used in studies of angiogenesis, mitogenesis of fibroblasts, neurite outgrowth in PC-12 cells, receptor binding, and tyrosine protein kinase phosphorylation.

EGF

Epidermal Growth Factor (EGF) has a profound effect on the differentiation of specific cells *in vivo* and is a potent mitogenic factor for a variety of cultured cells of both ectodermal and mesodermal origin.

Select your Gibco growth factor, cytokine, or chemokine at [thermofisher.com/growthfactors](https://www.thermofisher.com/growthfactors)

Cell culture reagents—antibiotics

Protecting your precious cultures or selecting for antibiotic resistance

Culture protection

Gibco™ cell culture reagents offer a wide range of antibiotics and antimycotics that can help control or eliminate cell culture contamination arising from bacteria, fungi, mycoplasma, and yeast.

Choose from these popular antibiotics:

Actinomycin D	Neomycin
Ampicillin	Penicillin
Gentamicin	Streptomycin (Pen Strep)
Kanamycin	Streptomycin

These and other antibiotics can be selected at [thermofisher.com/antibiotics](https://www.thermofisher.com/antibiotics)

Culture selection

Selection antibiotics for eukaryotic expression hosts provide:

- Fast selection of stable eukaryotic cell lines
- Selection in multiple organisms
- Dual selection



Eukaryotic selection antibiotics:

Blasticidin—resistance to blasticidin is conferred by the product of the *bsd* gene from *Aspergillus terreus*. Cell death occurs rapidly, and blasticidin-resistant stable mammalian cell lines can be generated in less than one week.

Geneticin (G-418)—Gibco™ Geneticin™ Selective Antibiotic exhibits higher purity so you can use less antibiotic and produce greater selection pressure with minimal toxicity from contaminants. Commonly known as G418 or G-418, the dominant-acting neomycin resistance gene used on common mammalian expression vectors contain elements derived from either transposons Tn601 (903) or Tn5. Stable colonies of mammalian cells expressing these resistance markers can be generated in 10 to 14 days.

Hygromycin B—a different mode of action from Geneticin Selective Antibiotic, Blasticidin S, or Zeocin™, makes Hygromycin B perfect for dual-selection experiments when used in conjunction with another selection agent. Resistance to Hygromycin B is conferred by the *E. coli* hygromycin resistance gene (*hyg* or *hph*).

Puromycin—resistance to puromycin is conferred by the puromycin N-acetyltransferase gene (*pac*). Puromycin has a fast mode of action, causing rapid cell death at low antibiotic concentrations. Puromycin-resistant stable mammalian cell lines can be generated in less than one week.

Zeocin—resistance to Zeocin is conferred by the *Sh ble* gene product. Effective in multiple cell types, eukaryotic expression vectors only need to carry one drug selection marker, which reduces the overall size of the vector and makes subcloning and transfection easier and more efficient.


Select your selection antibiotics at [thermofisher.com/selectionantibiotics](https://www.thermofisher.com/selectionantibiotics)

Cell culture reagents—cryopreservation

Upgrade home-brew freezing recipes to a complete, ready-to-use freezing medium

Cell lines in continuous culture are prone to genetic drift, finite cell lines are fated for senescence, all cell cultures are susceptible to microbial contamination, and even the best-run laboratories can experience equipment failure. Because an established cell line is a valuable resource and its replacement is expensive and time-consuming, it is vitally important that cells are frozen down and preserved for long-term storage.

Three ready-to-use Gibco freezing media are available for a range of different cell lines:



	Recovery™ Cell Culture Freezing Medium	Synth-a-Freeze™ Cryopreservation Medium	PSC Cryopreservation Kit
Description	Complete freezing medium for cryopreservation of a wide variety of mammalian cells	Intended for freezing and storing a variety of cell types	Cryopreservation medium and recovery supplement optimized to enable maximum pluripotent stem cell viability
Tested cell types	CHO-S, CHO-K1, HEK 293, Jurkat, NIH 3T3	Human keratinocytes, ESCs, MSCs, NSCs, other primary cell types	PSCs (iPSCs and ESC), PBMCs, iPSC-derived cardiomyocytes
Chemical composition	Contains FBS	Animal origin-free (AOF)	Xeno-free
Size	50 mL	50 mL	50 mL bottle + 5 mL vial

Select your cryopreservation reagent at [thermofisher.com/cryopreservation](https://www.thermofisher.com/cryopreservation)

Reach for recovery

Gibco Recovery Cell Culture Freezing Medium is a complete, ready-to-use freezing medium for a wide variety of commonly used mammalian cell lines

- No need to combine multiple products each time to make a home-brew freezing solution
- Improved viability of cells after thawing from cryostorage (Figure 5)




Figure 5. Comparative analysis of Recovery Cell Culture Freezing Medium vs. currently available cryopreservation technologies. For each of four different cell lines, the same number of cells was frozen in liquid nitrogen. Upon thawing, each cell sample was diluted into the same volume of a growth medium and seeded in 12-well plates (recovery was tested in eight different growth media). The initial seeding densities were 2.3×10^4 cells/well for adherent cell lines, and 6.0×10^3 cells/well for nonadherent cell lines, based on cell counts determined at the cryopreservation step. When the leading condition for a particular cell line had reached late log growth and none of the other conditions from that cell line had passed peak density, cells growing in all conditions for that cell line were harvested and the data analyzed. Thus, different cell lines were analyzed on different days: CHO, day 6; MDCK, day 6; Swiss 3T3, day 5; HeLa, day 6. Results were normalized to Recovery Cell Culture Freezing Medium and plotted with standard error bars ($n = 4$).

Mr. Frosty Freezing Container

Combine a ready-to-use Gibco cell freezing medium with a Thermo Scientific™ Mr. Frosty™ Freezing Container—optimal cryopreservation for superior cell recovery and viability post-thaw. Freeze cells in tubes from 1 to 5 mL using the Mr. Frosty Freezing Container. The simple-to-use system is designed to achieve a rate of cooling very close to $-1^\circ\text{C}/\text{minute}$, the optimal rate for cell preservation.

Features:

- Easy to use—requires only 100% isopropyl alcohol and a mechanical freezer
- Holds 12–18 tubes

Meet Mr. Frosty at thermofisher.com/mrfrosty

Reach for Recovery Freezing Medium at thermofisher.com/recovery

Cell culture reagents—matrices

The consistent base for your cell cultures

Gibco extracellular matrices, 3D scaffolds, and attachment proteins are essential tools for providing enhanced cell function resulting in more *in vivo*-like cell morphology, physiologically relevant environments, and better intercellular interactions.

- Proven technologies essential for cell attachment
- Lot-to-lot product consistency for ease of use and greater confidence
- Broad range of applications including liver and skin biology, stem cell research, neuroscience, and oncology

Geltrex™ Matrix—free of lactate dehydrogenase-elevating virus (LDEV), making Geltrex Matrix ideal for all types of cell culture and mouse *in vivo* cell culture research.

Collagen—one of the most widely used extracellular matrix proteins for cell culture, facilitating cell attachment, growth, differentiation, migration, and tissue morphogenesis.

Laminin—major glycoprotein component of basement membranes. Mouse laminin facilitates the attachment and expansion of embryonic stem cells (ESCs), induced pluripotent stem cells (iPSCs) and neural stem cells (NSCs).

Select these and other matrices at thermofisher.com/3D-cellculture




Cell culture—custom media

For times when cells require a unique media

Whether you need to customize an existing Gibco catalog product or have your own formulation that we can produce for you, we are committed to providing you with the products your research demands.

Not all projects are alike—each experiment can present unique needs and challenges. We are committed to providing you with cell culture products that are customized to your individual requirements:

- Tailor the format or package for special equipment
- Modify an existing formula (public or proprietary Gibco media)
- Derive your own formula from scratch
- Create media and packaging that perfectly fit your project



Learn more at [thermofisher.com/custommedia](https://www.thermofisher.com/custommedia)

Cell culture—equipment and instrumentation

Water baths

Precision redefined for your lab

Simplify workflows and maximize productivity with the new generation of Thermo Scientific™ Precision™ water baths. Enhanced designs and added features, combined with rugged construction and advanced microprocessor technology, make them the smart choice for your lab.

Select your water bath at [thermofisher.com/waterbath](https://www.thermofisher.com/waterbath)

Biological safety cabinets

Protection that never takes a day off

With Thermo Scientific™ biological safety cabinets, the certified performance and protection you get on Day 1 stays with you every day. With our proven reliability, ergonomics, and energy efficiency, the ideal choice is the one you can trust completely.

Select your biosafety cabinet at [thermofisher.com/biosafetycabinets](https://www.thermofisher.com/biosafetycabinets)

CO₂ incubators

Discovery thrives in a culture of confidence

More scientists worldwide trust their valuable cultures to Thermo Scientific™ CO₂ incubators than any other brand. They depend on proven reliability, outstanding contamination prevention, and optimal growth conditions. The reason is simple: our CO₂ incubators let you culture with confidence—day after day, year after year.

Select your CO₂ incubator at [thermofisher.com/co2](https://www.thermofisher.com/co2)

Centrifuges

Take your samples for a spin

Separations are a critical step in your workflows. That's why it's important to consider the centrifuge requirements and technical specifications for your applications, from selecting the appropriate speed and g-force to exploring the latest trends in centrifugation. Select your centrifuge at [thermofisher.com/centrifuges](https://www.thermofisher.com/centrifuges)

Freezers and refrigerators

Protect your samples

Proven cold-storage sample protection solutions from +4°C high-performance lab refrigerators to –196°C cryogenic freezers. Concentrate on your work without worrying about your valuable samples. Researchers worldwide protect more than two billion samples inside Thermo Scientific™ cold-storage equipment. Learn more about our temperature storage solutions at [thermofisher.com/cold](https://www.thermofisher.com/cold)

Automated cell counters

Accurate cell counts in just a few seconds

We offer two high-performance automated cell counters designed to meet the needs of any lab. The Invitrogen™ Countess™ II and Countess™ II FL Automated Cell Counters contain advanced autofocus and counting algorithms to allow you to quickly and accurately identify and count cells within a population. Applications include cell counting, viability, and transfection efficiency. Choose your Countess II instrument at [thermofisher.com/countessII](https://www.thermofisher.com/countessII)

Gibco Education

Culture your career

Gibco™ Education was developed out of our commitment to helping you achieve superior cell culture practices and outcomes every time you enter the lab. Our approach takes you from textbook concepts to successful cell culture practice in your lab.

Virtual labs

These courses offer a unique, interactive cell culture training experience. Learn everything from the fundamentals of cell culture to key techniques for primary and stem cell cultures and successful protein expression. As you progress through each lab, you can certify your knowledge by taking in-lab quizzes.

Explore the virtual labs to earn badges at thermofisher.com/gibcoeducation

Handbooks

The Gibco™ cell culture handbooks help your lab attain reproducible results every day.

Cell culture basics—download your copy at thermofisher.com/cellculturebasics

Protein expression—download your copy at thermofisher.com/expressionguide

Pluripotent stem cells—download your copy at thermofisher.com/pschandbook

Transfect your cells—download your copy at thermofisher.com/transfectionhandbook



Culture then transfect

Combine Lipofectamine transfection reagents with Gibco cell culture reagents for optimal efficiency and viability

Lipofectamine 3000

Single reagent kit with superior performance in the widest range of difficult- and easy-to-transfect cell types

- 10-fold higher transfection efficiency in difficult-to-transfect cells
- Lower cost per reaction vs. other transfection reagents

Learn more at [thermofisher.com/3000](https://www.thermofisher.com/3000)

Lipofectamine RNAiMAX

Superior transfection for your gene knockdown studies

- Transfect using less siRNA, resulting in lower cost
- Low cytotoxicity
- Works across multiple cell types

Learn more at [thermofisher.com/rnaimax](https://www.thermofisher.com/rnaimax)

Lipofectamine CRISPRMAX™

Achieve up to 85% cleavage efficiency with CRISPR-Cas9 protein delivery

- First optimized lipid nanoparticle transfection reagent for CRISPR-Cas9 protein delivery
- Gentle on cells
- Cost-effective

Learn more at [thermofisher.com/crisprmax](https://www.thermofisher.com/crisprmax)

Opti-MEM™ I

- Reduced-serum media ideal for use during cationic lipid transfections

Select your Opti-MEM media at [thermofisher.com/optimem](https://www.thermofisher.com/optimem)



Culture then express

Experience a fully optimized transient CHO expression system

The Gibco™ ExpiCHO™ Expression System combines Gibco cell culture media expertise with specialized cells and a unique transfection reagent to create an optimized system for protein expression.

- Achieve protein yields up to 3 g/L
- Express proteins in CHO cells from discovery to manufacture, providing an efficient and seamless drug development workflow
- Most cost-effective transient expression system available (more protein per dollar)

Discover the ExpiCHO Expression System at [thermofisher.com/expicho](https://www.thermofisher.com/expicho)




Figure 6. Recombinant protein titers in FreeStyle CHO, Expi293, and ExpiCHO systems. Expression levels of human IgG, rabbit IgG, and erythropoietin in the Gibco™ FreeStyle CHO, Gibco™ Expi293, and Gibco™ ExpiCHO transient expression systems are shown. Titrers in the ExpiCHO system range from 25x to 160x those of the FreeStyle CHO system, and 2x to 4x those obtained using the Expi293 system.

gibco



For more information, contact your local distributor



Austria: +43 (0)800 20 88 40 **Belgium:** +32 (0)056 260 260 **Denmark:** +45 70 27 99 20 **Finland:** +358 (0)9 8027 6280 **France:** +33 (0)3 88 67 14 14
Germany: +49 2304 932-5 **Netherlands:** +31 (0)20 487 70 00 **Norway:** +47 22 95 59 59 **Spain:** +34 902 239 303 **Sweden:** +46 31 352 32 00
UK: +44 (0)1509 555 500

ThermoFisher
SCIENTIFIC

For Research Use Only. Not for use in diagnostic procedures. © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Essential 8 is a trademark of Cellular Dynamics International, Inc. ScienCell is a trademark of ScienCell Research Laboratories. Lonza is a trademark of Lonza Group Ltd. B-27 is a trademark of Board of Trustees of Southern Illinois University. Zeocin is a trademark of InvivoGen. **COL31206 0416**