

WHY CHOOSE BIOCLEAN™ NITRILE ISOLATOR/RABS GLOVES?

Nitrile is a synthetic, non-solvent based, FDA approved polymer and is an ideal alternative to CSM, EPDM and latex, where the risk of latex allergies is a concern. With excellent anti-static properties, preventing the buildup of static electricity, nitrile is ideal for use with flammable liquids and powders. It can also be sanitised by Gamma Irradiation, Vapourised Hydrogen Peroxide (VHP) and Isopropyl Alcohol (IPA) and non-sterile options can also be washed, processed and packaged within a cleanroom environment, ensuring the gloves are an ultra-low contamination risk before being introduced into the isolator glove box.

ISOLATOR GLOVE MATERIAL COMPARISON

| Polymer | Sterilization Gamma | Repeated sterilization | | FDA | Cost | Comfort | Mechanical | Chemical |
|----------------------|------------------------|------------------------|-----|------------|----------|----------|------------|------------|
| | | Autoclave | VHP | Compliance | | | Properties | Resistance |
| Nitrile | | | | | | | | |
| CSM | | | | | | | | |
| EPDM | | | | | | O | | |
| EPDM Plus | | | | | | 0 | | |
| Natural Rubber Latex | | | | | | | | |
| Neoprene | | | | | | | | |
| | Poor | Fair | | Good | Very Goo | od | Excellent | |

YOUR CONFIDENCE, OUR PRIORITY

Ansell RABS and Isolator gloves have quality built in, because our quality control procedure is one of the most rigorous in the industry, including;








Final inspection controls;

Inspection for holes, visually and subsequently using a water and air pressure test









NITRILE PRODUCT PORTFOLIO OVERVIEW

| 1 | BioClean™ GGL | BioClean™ CGL | BioClean™ GGL30NITM9 | BioClean™ GSL | BioClean™ GSG10NIT80/ GSG10NIT85 | BioClean™ GSG10NITXLMA |
|----------------|----------------------------------|--------------------------------------|-------------------------|---------------|-------------------------------------|---------------------------|
| (BioClean | | 700 | | | | |
| Material | Nitrile | Nitrile | Nitrile | Nitrile | Nitrile | e/ PCP |
| Style | Glove | Glove | Mitten | Sleeve | Sleeve Glo | ve System |
| Sterility | Sterile & Cleanroom Laundered | Non-Sterile & Cleanroom Laundered | Sterile or Non-Sterile | Sterile | Ster | rile |
| Cuff Thickness | 0.50mm/20mil | 0.50mm/20mil | 0.50mm/20mil | 0.50mm/20mil | 0.50mm/20mil | 0.50mm/20mil |
| Port Size | 8, 10, 12" | 10, 12" | 10, 12" | 6, 8, 10, 12" | 8, 1 | 10" |
| Length | 840mm/33" | 840mm/33" | 840mm/33" | 660mm/26" | 900mm / 35" | 995 mm / 39" |
| Hand Size | 9.75/L | 9.75/L | 9.75/L | 9.75/L | 8/M, 8 | 3.5/M |

NITRILE MATERIAL BENEFITS

- FDA approved polymer
- Excellent anti-static properties
- Performs well with VHP or IPA
- · Superior comfort and dexterity
- Excellent chemical resistance

CHEMICAL PERMEATION RESULTS TABLE - BIOCLEAN™ GGL, CGL, GHG & CHG

*ASTM F739 - Breakthrough of the test chemical is deemed to have occurred when the permeation rate has reached $0.1 \,\mu\text{g/cm}^2$ /min. **EN 16523-1: 2015 (formerly EN374-3) - Breakthrough of the test chemical is deemed to have occurred when the permeation rate has reached $1.0 \,\mu\text{g/cm}^2$ /min.

| PERMEATION BREAKTHROUGH TIMES | 0 | 1 | 2 | 3 | 4 | 5 | 6 |
|----------------------------------|-----------------|-------------------|-------|-------------------|---------|-----------------|------|
| | <10 | 10-30 | 30-60 | 60-120 | 120-240 | 240-480 | >480 |
| (MINUTES) | Not Recommended | Splash Protection | | Medium Protection | | High Protection | |

| | 0, | | STANDARD | | |
|---|-----|-----------|----------|---------|--|
| CHEMICAL NAME | | CAS | ASTM* | EN** | |
| 1-BUTANOL (BUTYL ALCOHOL) | 100 | 71-36-3 | >480 | >480 | |
| 1-PROPANOL (PROPYLALCOHOL, N-PROPANOL) | 100 | 71-23-8 | 240-480 | 240-480 | |
| ACETIC ACID, GLACIAL † | 99 | 64-19-7 | 110 | 107 | |
| ACETONE (2-PROPANONE) | 100 | 67-64-1 | 7 | 7 | |
| CITRIC ACID | 100 | 77-92-9 | >480 | >480 | |
| CYCLOHEXANE | 100 | 110-82-7 | 240-480 | 240-480 | |
| ETHANOL (ETHYLALCOHOL) | 100 | 64-17-5 | 278 | 364 | |
| FORMALDEHYDE † | 37 | 50-00-0 | >480 | >480 | |
| HEPTANE (N-HEPTANE) † | 100 | 142-82-5 | 240-480 | >480 | |
| HEXANE (N-HEXANE) | 100 | 110-54-3 | >480 | >480 | |
| HYDROCHLORIC ACID | 37 | 7647-01-0 | >480 | >480 | |
| HYDROGEN PEROXIDE | 37 | 7722-84-1 | >480 | >480 | |
| ISOBUTANOL (ISOBUTYLALCOHOL) | 100 | 78-83-1 | 240-480 | 240-480 | |
| ISOPROPYLALCOHOL (IPA, ISOPROPANOL, 2-PROPANOL) | 100 | 67-63-0 | >480 | >480 | |
| METHANOL † | 100 | 67-56-1 | 56 | 57 | |
| METHYL ETHYL KETONE (2-BUTANONE, MEK) | 100 | 78-93-3 | <10 | <10 | |
| PHOSPHORIC ACID | 45 | 7664-38-2 | >480 | >480 | |
| SODIUM HYDROXIDE (NAOH) † | 50 | 1310-73-2 | >480 | >480 | |
| SODIUM HYPOCHLORITE | 8.5 | 7681-52-9 | >480 | >480 | |
| SULFURIC ACID | 50 | 7664-93-9 | >480 | >480 | |
| SULFURIC ACID | 96 | 7664-93-9 | <10 | <10 | |

CHEMICAL PERMEATION RESULTS TABLE - BIOCLEAN™ GGL, CGL, GHG & CHG

*ASTM F739 - Breakthrough of the test chemical is deemed to have occurred when the permeation rate has reached 0.1 µg/cm² /min.

**EN 16523-1: 2015 (formerly EN374-3) - Breakthrough of the test chemical is deemed to have occurred when the permeation rate has reached
1.0 µg/cm² /min.

| THERE BARTY DISINIFICANT REANING | STANDARD | | |
|----------------------------------|----------|---------|--|
| THIRD PARTY DISINFECTANT BRANDS | ASTM* | EN** | |
| DECON-CLEAN® | 240-480 | 240-480 | |
| DECON-SPORE 200® PLUS | 240-480 | 240-480 | |
| KLERCIDE™ CR BIOCIDE S | 240-480 | 240-480 | |
| KLERCIDE™ Y | 240-480 | 240-480 | |
| LPH® SE | 240-480 | 240-480 | |
| SPORE-KLENZ® | 240-480 | 240-480 | |
| VESPHENE® IISE | 240-480 | 240-480 | |

[†]Certified test result

Please see product validation pack for full permeation results. When a number is listed in a cell, this means that an actual test has been performed. The number is showing the permeation time in minutes. When a cell is coloured, with no number, the permeation time is based on extrapolation issued from AnsellGUARDIAN®. When a cell is white, no data is available.

The chemical permeation results table is related to the barrier performance of certain personal protective equipment (PPE) against the chemicals. This information is intended to enable the Health and Safety professional at your organisation make more informed decisions about the Ansell PPE that may offer the greatest protection in the intended circumstances and assist with carrying out a risk assessment for your organisation. We wish to highlight that permeation times do not equate to safe wear time. Safe wear time may vary depending on whether the PPE is donned correctly, the surrounding temperature, the chemicals' toxicity, and other factors. Permeation information offered here is limited to the main protective material. Permeation times may vary around seams, zips, visors or any other joins or components of the PPE. It is the responsibility of your organisation's Health and Safety professional to undertake a risk assessment before choosing the appropriate PPE for the task at hand. If you want to discuss any aspect in detail, please contact us.

Estimations of the barrier properties of PPE are based on currently available data and extrapolations from laboratory test results and information regarding the chemicals' composition. Synergistic effects of mixing chemicals have not been accounted for. Estimations are subject to change if new testing is carried out or new information is available providing better grounds for extrapolations. For these reasons, any information in this report is provided for informational purposes only and Ansell fully disclaims any liability including warranties related to any statement contained herein.

