


Product Description :	Face mask PP, white, 2 plies, with ear loop
Product Code :	15451932 (17-100-896)
Raw Material :	18g/m ² + 21g/m ² Polypropylene
Size :	17.5 x 9.5 cm
Packaging :	50 pcs/bag, 10 bags/box
Carton Dimension :	37 x 22 x 20 cm



Item	Length	Allowance	Remark
A	17.5cm	+/- 1cm	
B	9.5cm	+/- 1cm	
C	11.5cm	+/- 1cm	
D	16cm	+/- 1cm	

© 2017 Thermo Fisher Scientific Inc. All rights reserved.
 Trademarks used are owned as indicated at fishersci.com/trademarks.

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
 Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500